

Avant™ buffer pouches

Convenient, pure, fast buffer preparation

Avansta's **Avant buffer pouches** are perfect for standardizing your electrophoresis and Western blotting applications. Dissolve a pouch of pre-measured molecular biology-grade chemicals in 500 ml of deionized water to obtain a ready-to-use 1x buffer without any further preparation. Save precious time in the laboratory and generate reproducible results with Avant buffers.



Advantages

- **QUALITY** – molecular biology-grade buffers of highest purity
- **CONVENIENT** – pre-measured powder in sealed pouches
- **FAST** – dissolve and use right away
- **STANDARDIZED** – guaranteed reproducible results for your electrophoresis gels and Western blots

Ordering Information

Catalog Number	Product	Size
R-01038-020	PBS (concentration 1x at 500 ml)	20 pouches
R-01039-020	TBS (concentration 1x at 500 ml)	20 pouches
R-01037-020	Tris-Glycine PAGE running buffer (concentration 1x at 500 ml)	20 pouches
R-01036-020	Tris-Glycine SDS PAGE running buffer (concentration 1x at 500 ml)	20 pouches

Related Products

K-12042-D10	WesternBright® Quantum® Western Blotting HRP Substrate	For 1000 cm ² membrane	100 ml
K-12045-D20	WesternBright® ECL Western Blotting HRP Substrate	For 2000 cm ² membrane	200 ml

Avansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@avansta.com

Product information: www.avansta.com/products/Avant_pouches

Copyright © 2010-2020 Avansta. All rights reserved. The Avansta logo, Quantum® and WesternBright® are registered trademarks of the Company. Avant™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

www.avansta.com

