


Blot Development Folders

Transparent plastic supports for blot handling and imaging

Avansta's **Blot Development Folders** are transparent plastic folders that hold your blot while it is exposed to film or CCD. These folders keep your blot flat for perfect film exposures, avoiding the wrinkling or bunching that can occur with plastic wrap.

Advantages


- **CLEAR IMAGES** – avoid the wrinkling or bunching that can happen with plastic wrap
- **SUITABLE FOR ALL BLOTS** – development folders are suitable for chemiluminescence and fluorescence imaging
- **SAFE BLOT HANDLING** – move your blot without scratching or putting pressure on the surface
- **CONVENIENT SIZE** – image up to 2 blots simultaneously in one folder




Avansta Solutions

Want the cleanest possible images of your HRP-imaged blot?

Use Avansta's Blot Development Folders, which keep your blot neat and flat during film exposures and help prevent touch artifacts.



Blot Development Folders provide a clean, easy alternative to plastic wrap.



Ordering Information

Catalog Number	Product		Size
L-07020-025	Blot Development Folders	5x7 inches (12.7x17.8 cm)	25/pack
L-07020-100	Blot Development Folders	5x7 inches (12.7x17.8 cm)	100/pack

Related Products

K-12045-D20	WesternBright® ECL Western Blotting HRP Substrate	For 2000 cm ² membrane	200 ml
L-07014-100	LucentBlue™ X-ray Film	8x10 inches (20.3x25.4 cm)	100 sheets
L-07019-001	Film Cassette	8x10 inches (20.3x25.4 cm)	1 each
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm	8x10 cm	10 sheets

Avansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@avansta.com

Product information: www.avansta.com/products/Development_Folders

Copyright © 2014-2020 Avansta. All rights reserved. The Avansta logo and WesternBright® are registered trademarks of the Company. LucentBlue™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

www.avansta.com

