

Blotting Sponge Pads

More reproducible gel-to-membrane electrophoretic transfers

Advansta's **Blotting Sponge Pads** maintain secure contact between gel and membrane in transfer sandwiches and cassettes for consistent transfers. The Blotting Sponge Pads resist compression and are available in two sizes to work with multiple transfer setups and membrane sizes.



Advantages

- **RESIST COMPRESSION** – maintain proper pressure in transfer sandwiches without adding extra blotting papers
- **RETAIN THICKNESS** – can be reused over multiple experiments
- **CONVENIENT** – two sizes compatible with a variety of transfer apparatuses
- **STANDARDIZED** – uniform flat construction ensures reproducible and consistent transfers

Ordering Information

Catalog Number	Product		Size
L-07129-004	Blotting Sponge Pads	8 x 11 cm	4 sponges
L-07130-004	Blotting Sponge Pads	9 x 15 cm	4 sponges

Related Products

K-12042-D10	WesternBright® Quantum® Western Blotting HRP Substrate	For 1000 cm ² membrane	100 ml
L-07045-060	Pre-cut Blotting Paper	7 x 9 cm	60 sheets
L-07046-060	Pre-cut Blotting Paper	8 x 10 cm	60 sheets
L-08002-010	Pre-cut WesternBright® NC, 0.45 µm	8 x 10 cm	10 sheets
L-08003-010	Pre-cut WesternBright® NC, 0.22 µm	8 x 10 cm	10 sheets
L-08001-010	Pre-cut WesternBright® PVDF-FL	7 x 9 cm	10 sheets
L-08012-010	Pre-cut WesternBright® PVDF-FL	10 x 15 cm	10 sheets
L-08004-010	Pre-cut WesternBright® PVDF-CL	7 x 9 cm	10 sheets
L-08011-010	Pre-cut WesternBright® PVDF-CL	10 x 15 cm	10 sheets
R-03090-D25	FLASHBlot™ Transfer Buffer, 50x concentrate	Sufficient for up to 50 blots	250 ml
R-03090-D50	FLASHBlot™ Transfer Buffer, 50x concentrate	Sufficient for up to 100 blots	500 ml

Advansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com

Product information: www.advansta.com/products/Sponge-Pads

Copyright © 2015-2020 Advansta. All rights reserved. The Advansta logo, Quantum® and WesternBright® are registered trademarks of the Company. FLASHBlot™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

www.advansta.com

