


# LucentBlue™ X-ray film

High sensitivity X-ray film for chemiluminescent Western blots

**LucentBlue** is a sensitive film ideal for imaging chemiluminescent Western blots and is guaranteed to perform with WesternBright® ECL, WesternBright Quantum® and WesternBright Sirius®. LucentBlue is compatible with horseradish peroxidase and alkaline phosphatase detection reagents, as well as autoradiography. In addition to Western blots, LucentBlue can be used for Northern and Southern blot detection, and other techniques requiring film.



## Advantages

- **PERFORMANCE** – excellent images when used with WesternBright ECL, WesternBright Quantum, and WesternBright Sirius
- **VALUE** – sensitivity equal to more expensive films
- **FLEXIBILITY** – detect radioisotopes in addition to chemiluminescence

## Ordering Information

Catalog Number	Product		Size
L-07014-100	LucentBlue™ X-ray film	8 x 10" sheets	100 sheets
L-07013-100	LucentBlue™ X-ray film	5 x 7" sheets	100 sheets

## Related Products

K-12045-D20	WesternBright® ECL Western Blotting HRP Substrate	For 2000 cm <sup>2</sup> membrane	200 ml
K-12042-D20	WesternBright® Quantum® Western Blotting HRP Substrate	For 2000 cm <sup>2</sup> membrane	200 ml
K-12043-D20	WesternBright® Sirius® Western Blotting HRP Substrate	For 2000 cm <sup>2</sup> membrane	200 ml
R-05071-500	Goat-anti-mouse HRP-conjugated secondary antibody		500 µl
R-05072-500	Goat-anti-rabbit HRP-conjugated secondary antibody		500 µl
L-08001-010	Low-Fluorescence PVDF Transfer Membrane	7 x 9 cm	10 sheets
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm	8 x 10 cm	10 sheets
L-08003-010	Nitrocellulose Transfer Membrane 0.22 µm	8 x 10 cm	10 sheets

## Advanta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: [sales@advanta.com](mailto:sales@advanta.com)

Product information: [www.advanta.com/products/lucentblue](http://www.advanta.com/products/lucentblue)

Copyright © 2010-2020 Advanta. All rights reserved. The Advanta logo, Quantum®, Sirius® and WesternBright® are registered trademarks of Advanta Corporation. LucentBlue™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

[www.advanta.com](http://www.advanta.com)

