

# SpectraDye™ Secondary Antibodies

## Fluorescently labeled antibodies

Advansta offers a broad selection of **SpectraDye Secondary Antibodies** with a variety of conjugate types, including visible and near infrared (NIR) fluorescent dyes. Choose the label and species best suited for your application. Each high quality antibody offers unparalleled sensitivity and performance. Additionally, achieve fast and convenient workflow when using the dyes in conjunction with Advansta's WesternBright® MCF system for multi-color fluorescent Westerns.



## Advantages

- **PERFORMANCE** – strong signal and no cross-reactivity
- **CONVENIENCE** – each antibody is demonstrated to provide optimal results with Advansta's Western blotting systems
- **FLEXIBILITY** – choose the label and target species best suited to your application

	Excitation maximum (nm)	Emission maximum (nm)
■ Dye-490	491	515
■ Dye-550	551	565
■ Dye-650	653	672
■ Dye-IR700	690	709
■ Dye-IR800	783	800

## Ordering Information

Product	Fluorescent label				
	490	550	650	IR700	IR800
Goat-anti-rabbit	R-05759-250	R-05760-250	R-05761-250	R-05054-250	R-05060-250
Goat-anti-mouse	R-05762-250	R-05763-250	R-05764-250	R-05055-250	R-05061-250
Goat-anti-human	R-05765-250	R-05766-250	R-05767-250	R-05056-250	R-05062-250
Goat-anti-rat	R-05768-250	R-05769-250	R-05770-250	R-05058-250	R-05064-250
Goat-anti-chicken	R-05771-250	R-05772-250	R-05773-250	R-05057-250	R-05063-250
Goat-anti-guinea pig	R-05774-250	R-05775-250	R-05776-250	R-05059-250	R-05065-250
Donkey-anti-goat	R-05777-250	R-05778-250	R-05779-250	R-05780-250	R-05781-250
Donkey-anti-mouse	R-05789-250	R-05790-250	R-05791-250	R-05783-250	R-05784-250
Donkey-anti-rabbit	R-05792-250	R-05793-250	R-05794-250	R-05785-250	R-05786-250

## Related Products

R-05051-250	Goat-anti-rabbit IgG APC
R-05052-250	Goat-anti-mouse IgG RPE

\*Size 250 µl for all products listed above.

## Advansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: [sales@advansta.com](mailto:sales@advansta.com)

Product information: [www.advansta.com/Fluorescent\\_secondary\\_antibodies](http://www.advansta.com/Fluorescent_secondary_antibodies)

Copyright © 2010–2020 Advansta. All rights reserved. The Advansta logo and WesternBright® are registered trademarks of Advansta Corporation. SpectraDye™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

[www.advansta.com](http://www.advansta.com)

